SAFE PERCUTANEOUS D R E N T E DRAINAGE KIT UN

For all-purposes drainage procedures

13 60

10.00



Η

MULTI

J SELDINGER

FORTY



PULL & PUSH

Thanks to double integrated valve drainage procedure can be completed by mean of simple "pull and push" manoeuvres. Three-way stopcock no longer needed.

PULL

PUSH

EXTENSION LINE

By connecting an extension line it improves patient comfort in case of extended use of the device and more flexibility for users.

DOUBLE INTEGRATED ONE-WAY VALVE

W FF. ET. SF. IF. OF. S.

designed to improve comfort and safety of your drainage procedures. The valves ensure patient and users safety by reducing risk of external environment contamination, pneumothorax and fluid leakage.

Veress needle with luer lock port for primary positioning check.

VERRES NEEDLE

retractable tip to achieve safe and efficient access. The blunt tip retracts during penetration and extends when the needle enters free space.

POLYURETHANE CATHETER:

totally radiopaque, anti kinking and cut-resistant. It is biocompatible and thermosensitive which makes it the first choice for long-term drainage. Larger inner diameter increases drainage flow rate. The best for pneumothorax treatment (primary and secondary), malignant or chronic pleural effusions and uncomplicated empyema.





PCS/BOX

STANDARD

Percutaneous drainage Kit - CH/Fr 9	10800	6
Percutaneous drainage Kit - CH/Fr 12	10801	6
Percutaneous drainage Kit with extension line - CH/Fr 9	10802	6
Percutaneous drainage Kit with extension line - CH/Fr 12	10803	6
Percutaneous drainage Kit with extension line - CH/Fr 15	10830	6

EXTRA LONG*

Percutaneous drainage Kit with extension line XL- CH/Fr 9	10820	
Percutaneous drainage Kit with extension line XL- CH/Fr 12	10821	6

PCS/BOX

*50% Longer than standard


Particularly suitable for abscess, nephrostomy pericardiocentesis and miscellaneous drainage procedures.



M			
	U		

Percutaneous drainage Kit "MULTI" CH/Fr 9	10816	6
Percutaneous drainage Kit "MULTI" CH/Fr 12	10817	6
Percutaneous drainage Kit "MULTI" with extension line - CH/Fr 9	10818	6
Percutaneous drainage Kit "MULTI" with extension line - CH/Fr 12	10819	6
Percutaneous drainage Kit "MULTI" with extension line - CH/Fr 15	10831	6

CODE PCS/BOX



91.

6. 8. 6.

Pig tail provided with five inner drainage holes, to prevent irritation of surrounding tissues and occlusion of the catheter The most precise positioning through a fast-track Seldinger technique, all together with the unmatched UNICO's features.



J SELDINGER

UNICO-J Percutaneous Drainage Seldinger - CH/Fr 9	10823	6
UNICO-J Percutaneous Drainage Seldinger - CH/Fr 12	10824	6

CODE PCS/BOX

Integrated

dilator

Verres needle

Atraumatic

J-type guide wire



ONE SHOT PUNCTURE SET

For pleural drainage of liquids and air, particularly suitable for: pneumothorax, effusion and low viscosity empyema

The forty centimeters percutaneous drainage design for apical positioning to ease liquids and air evacuation

FORTYCODEPCS/BOXPercutaneous drainage Kit with extension line - CH/Fr 9108106Percutaneous drainage Kit with extension line - CH/Fr 12108126

BASE

Unico Base - gauge 13/ 2,5 mm	10815	20

CODE

PCS/BOX

ACCESSORIES

	GODE	P03/ B0A
Drainage bag	12019	25
Catheter fixing kit	10546	50
M connector for chest drain	10548	50
F connector for chest drain	10504	20

10546 - Catheter

10548 - M connector for chest drain

fixing kit

INTEGRATED ONE-WAY VALVE to prevent liquid back flow

BUILT-IN Venting valve

10504 - F connector for chest drain

REDAX™

REDAX - Via Galileo Galilei, 18 46025 Poggio Rusco (MN) ITALY Tel. +39 0386 830582 Fax +39 0386 51898 E-mail: info@redax.it - www.redax.it



UNICO